

Protocol: adult\_upperExtremities\_Wrist-Supine Feet 1st

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>Gradient Echo</i>
Coil Configuration	<i>GEM Flex Small Array</i>	Imaging Options	<i>Seq, Fast</i>
Plane	<i>3-PLANE</i>	<b>SCAN RANGE</b>	
Series Description	<i>3-Plane Loc</i>	FOV	<i>36.0</i>
<b>SCAN TIMING</b>		Slice Thickness	<i>5.0</i>
Number of Echoes	<i>1</i>	Slice Spacing	<i>5.0</i>
<b>IMAGE ENHANCE</b>		<b>ACQ TIMING</b>	
Filter Choice	<i>None</i>	Freq	<i>256</i>
<b>GATING/TRIGGER</b>		Phase	<i>128</i>
Auto Trigger Type	<i>Off</i>	Freq DIR	<i>Unswap</i>
<b>MULTI-PHASE</b>		NEX	<i>1.00</i>
Seperate Series	<i>0</i>	Phase FOV	<i>1.00</i>
Mask Phase	<i>0</i>	Auto Shim	<i>Auto</i>
Mask Pause	<i>0</i>	Phase Correction	<i>No</i>
<b>DIFFUSION</b>		<b>FMRI</b>	
Recon All Images	<i>On</i>	PSD Trigger	<i>Internal</i>
<b>CONTRAST</b>		View Order	<i>Bottom/Up</i>
Contrast Yes/No	<i>No</i>	# of Repetitions REST	<i>0</i>
		# of Repetitions ACTIVE	<i>0</i>
		<b>SAT</b>	
		Tag Type	<i>None</i>
		<b>TRICKS</b>	
		Pause On/Off	<i>On</i>
		Auto Subtract	<i>0</i>
		Auto SCIC	<i>Off</i>

3-Plane Loc

3-Plane Loc

Protocol: adult\_upperExtremities\_Wrist-Supine Feet 1st

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>Gradient Echo</i>
Coil Configuration	<i>GEM Flex Small Array</i>	Imaging Options	<i>Seq, Fast</i>
Plane	<i>3-PLANE</i>	<b>SCAN RANGE</b>	
Series Description	<i>3-Plane Loc</i>	FOV	<i>15.0</i>
<b>SCAN TIMING</b>		Slice Thickness	<i>5.0</i>
Number of Echoes	<i>1</i>	Slice Spacing	<i>5.0</i>
<b>IMAGE ENHANCE</b>		<b>ACQ TIMING</b>	
Filter Choice	<i>None</i>	Freq	<i>256</i>
<b>GATING/TRIGGER</b>		Phase	<i>128</i>
Auto Trigger Type	<i>Off</i>	Freq DIR	<i>Unswap</i>
<b>MULTI-PHASE</b>		NEX	<i>1.00</i>
Seperate Series	<i>0</i>	Phase FOV	<i>1.00</i>
Mask Phase	<i>0</i>	Auto Shim	<i>Auto</i>
Mask Pause	<i>0</i>	Phase Correction	<i>No</i>
<b>DIFFUSION</b>		<b>FMRI</b>	
Recon All Images	<i>On</i>	PSD Trigger	<i>Internal</i>
<b>CONTRAST</b>		View Order	<i>Bottom/Up</i>
Contrast Yes/No	<i>No</i>	# of Repetitions REST	<i>0</i>
		# of Repetitions ACTIVE	<i>0</i>
		<b>SAT</b>	
		Tag Type	<i>None</i>
		<b>TRICKS</b>	
		Pause On/Off	<i>On</i>
		Auto Subtract	<i>0</i>
		Auto SCIC	<i>Off</i>

3-Plane Loc

3-Plane Loc

Protocol: adult\_upperExtremities\_Wrist-Supine Feet 1st

Ax T1	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>		Ax T1
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>	
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>FSE-XL</i>	
	Coil Configuration	<i>GEM Flex Small Array</i>	Imaging Options	<i>EDR, TRF, Fast</i>	
	Plane	<i>OBLIQUE</i>	<b>SCAN RANGE</b>		
	Series Description	<i>Ax T1</i>	FOV	<i>10.0</i>	
	<b>SCAN TIMING</b>		Slice Thickness	<i>3.0</i>	
	Flip Angle	<i>120</i>	Slice Spacing	<i>0.6</i>	
	TE	<i>Min Full</i>	<b>ACQ TIMING</b>		
	Number of Echoes	<i>1</i>	Freq	<i>320</i>	
	TR	<i>534.0</i>	Phase	<i>192</i>	
	Echo Train Length	<i>3</i>	Freq DIR	<i>R/L</i>	
	Receiver Bandwidth	<i>31.25</i>	Fat Shift DIR	<i>Normal (L)</i>	
	<b>IMAGE ENHANCE</b>		NEX	<i>3.00</i>	
	Filter Choice	<i>None</i>	# of Acq. Before Pause	<i>0</i>	
	<b>GATING/TRIGGER</b>		Phase FOV	<i>1.00</i>	
	Auto Trigger Type	<i>Off</i>	Auto Shim	<i>Auto</i>	
	<b>FMRI</b>		Phase Correction	<i>Yes</i>	
	PSD Trigger	<i>Internal</i>	<b>USER CVS</b>		
	View Order	<i>Bottom/Up</i>	User CV8	<i>2.00</i>	
	# of Repetitions REST	<i>0</i>	User CV9	<i>1.00</i>	
	# of Repetitions ACTIVE	<i>0</i>	User CV22	<i>1.00</i>	
	<b>SAT</b>		TR Min	<i>350.0</i>	
	SAT Location	<i>S</i>	TR Max	<i>800.0</i>	
	Tag Type	<i>None</i>	<b>MULTI-PHASE</b>		
<b>TRICKS</b>		Seperate Series	<i>0</i>		
Pause On/Off	<i>On</i>	Mask Phase	<i>0</i>		
Auto Subtract	<i>0</i>	Mask Pause	<i>0</i>		
Auto SCIC	<i>Off</i>	<b>DIFFUSION</b>			
		Recon All Images	<i>On</i>		
		<b>CONTRAST</b>			
		Contrast Yes/No	<i>No</i>		

Protocol: adult\_upperExtremities\_Wrist-Supine Feet 1st

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Feet First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	FSE-XL
Coil Configuration	GEM Flex Small Array	Imaging Options	Seq, EDR, TRF, Fast, FR, IrP
Plane	OBLIQUE		
Series Description	Ax STIR		
SCAN TIMING		SCAN RANGE	
Flip Angle	160	FOV	10.0
TE	40.0	Slice Thickness	3.0
Number of Echoes	1	Slice Spacing	0.6
TR	4353.0		
TI	140	ACQ TIMING	
Echo Train Length	8	Freq	288
Receiver Bandwidth	35.71	Phase	192
		Freq DIR	R/L
		Fat Shift DIR	Normal (L)
		NEX	3.00
		Phase FOV	1.00
		Auto Shim	Auto
		Phase Correction	Yes
IMAGE ENHANCE		USER CVS	
Filter Choice	None	User CV8	1.00
		User CV22	1.00
		TR Min	2500.0
		TR Max	8000.0
GATING/TRIGGER		MULTI-PHASE	
Auto Trigger Type	Off	Seperate Series	0
		Mask Phase	0
		Mask Pause	0
FMRI		DIFFUSION	
PSD Trigger	Internal	Recon All Images	On
View Order	Bottom/Up		
# of Repetitions REST	0	CONTRAST	
# of Repetitions ACTIVE	0	Contrast Yes/No	No
SAT			
SAT Location	S		
Tag Type	None		
TRICKS			
Pause On/Off	On		
Auto Subtract	0		
Auto SCIC	Off		

Ax STIR

Ax STIR

Protocol: adult\_upperExtremities\_Wrist-Supine Feet 1st

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>FSE-XL</i>
Coil Configuration	<i>GEM Flex Small Array</i>	Imaging Options	<i>NPW, EDR, TRF, Fast</i>
Plane	<i>OBLIQUE</i>	<b>SCAN RANGE</b>	
Series Description	<i>Cor T1</i>	FOV	<i>10.0</i>
<b>SCAN TIMING</b>		Slice Thickness	<i>3.0</i>
Flip Angle	<i>120</i>	Slice Spacing	<i>0.0</i>
TE	<i>Min Full</i>	<b>ACQ TIMING</b>	
Number of Echoes	<i>1</i>	Freq	<i>288</i>
TR	<i>581.0</i>	Phase	<i>224</i>
Echo Train Length	<i>3</i>	Freq DIR	<i>S/I</i>
Receiver Bandwidth	<i>35.71</i>	Fat Shift DIR	<i>Normal (S)</i>
<b>IMAGE ENHANCE</b>		NEX	<i>3.00</i>
Filter Choice	<i>None</i>	Auto Shim	<i>Auto</i>
<b>GATING/TRIGGER</b>		Phase Correction	<i>Yes</i>
Auto Trigger Type	<i>Off</i>	<b>USER CVS</b>	
<b>FMRI</b>		User CV8	<i>2.00</i>
PSD Trigger	<i>Internal</i>	User CV9	<i>1.00</i>
View Order	<i>Bottom/Up</i>	User CV22	<i>1.00</i>
# of Repetitions REST	<i>0</i>	TR Min	<i>350.0</i>
# of Repetitions ACTIVE	<i>0</i>	TR Max	<i>800.0</i>
<b>SAT</b>		<b>MULTI-PHASE</b>	
Tag Type	<i>None</i>	Seperate Series	<i>0</i>
<b>TRICKS</b>		Mask Phase	<i>0</i>
Pause On/Off	<i>On</i>	Mask Pause	<i>0</i>
Auto Subtract	<i>0</i>	<b>DIFFUSION</b>	
Auto SCIC	<i>Off</i>	Recon All Images	<i>On</i>
		<b>CONTRAST</b>	
		Contrast Yes/No	<i>No</i>

Cor T1

Cor T1

Protocol: adult\_upperExtremities\_Wrist-Supine Feet 1st

Cor T2W	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>		Cor T2W
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>	
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>FSE-XL</i>	
	Coil Configuration	<i>GEM Flex Small Array</i>	Imaging Options	<i>NPW, EDR, TRF, Fast, FR</i>	
	Plane	<i>OBLIQUE</i>	<b>SCAN RANGE</b>		
	Series Description	<i>Cor T2W</i>	FOV	<i>10.0</i>	
	<b>SCAN TIMING</b>		Slice Thickness	<i>3.0</i>	
	Flip Angle	<i>160</i>	Slice Spacing	<i>0.0</i>	
	TE	<i>77.0</i>	<b>ACQ TIMING</b>		
	Number of Echoes	<i>1</i>	Freq	<i>288</i>	
	TR	<i>2588.0</i>	Phase	<i>224</i>	
	Echo Train Length	<i>16</i>	Freq DIR	<i>S/I</i>	
	Receiver Bandwidth	<i>35.71</i>	Fat Shift DIR	<i>Normal (S)</i>	
	<b>IMAGE ENHANCE</b>		NEX	<i>4.00</i>	
	Filter Choice	<i>None</i>	Auto Shim	<i>Auto</i>	
	<b>GATING/TRIGGER</b>		Phase Correction	<i>Yes</i>	
	Auto Trigger Type	<i>Off</i>	<b>USER CVS</b>		
	<b>FMRI</b>		User CV8	<i>2.00</i>	
	PSD Trigger	<i>Internal</i>	User CV9	<i>1.00</i>	
	View Order	<i>Bottom/Up</i>	User CV22	<i>1.00</i>	
	# of Repetitions REST	<i>0</i>	TR Min	<i>2500.0</i>	
	# of Repetitions ACTIVE	<i>0</i>	TR Max	<i>6000.0</i>	
	<b>SAT</b>		<b>MULTI-PHASE</b>		
	Tag Type	<i>None</i>	Seperate Series	<i>0</i>	
	<b>TRICKS</b>		Mask Phase	<i>0</i>	
Pause On/Off	<i>On</i>	Mask Pause	<i>0</i>		
Auto Subtract	<i>0</i>	<b>DIFFUSION</b>			
Auto SCIC	<i>Off</i>	Recon All Images	<i>On</i>		
		<b>CONTRAST</b>			
		Contrast Yes/No	<i>No</i>		

Protocol: adult\_upperExtremities\_Wrist-Supine Feet 1st

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Feet First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	FSE-XL
Coil Configuration	GEM Flex Small Array	Imaging Options	NPW, Seq, EDR, TRF, Fast, ZIP512, FR, IrP
Plane	OBLIQUE		
Series Description	Cor STIR		
SCAN TIMING		SCAN RANGE	
Flip Angle	160	FOV	10.0
TE	34.0	Slice Thickness	3.0
Number of Echoes	1	Slice Spacing	0.0
TR	3551.0		
TI	130	ACQ TIMING	
Echo Train Length	10	Freq	256
Receiver Bandwidth	35.71	Phase	160
		Freq DIR	S/I
		Fat Shift DIR	Normal (S)
IMAGE ENHANCE		NEX	4.00
Filter Choice	None	Auto Shim	Auto
		Phase Correction	Yes
GATING/TRIGGER		USER CVS	
Auto Trigger Type	Off	User CV8	1.00
		User CV22	1.00
FMRI		TR Min	2500.0
PSD Trigger	Internal	TR Max	8000.0
View Order	Bottom/Up		
# of Repetitions REST	0	MULTI-PHASE	
# of Repetitions ACTIVE	0	Seperate Series	0
		Mask Phase	0
SAT		Mask Pause	0
Tag Type	None	DIFFUSION	
		Recon All Images	On
TRICKS		CONTRAST	
Pause On/Off	On	Contrast Yes/No	No
Auto Subtract	0		
Auto SCIC	Off		

Cor STIR

Cor STIR

Protocol: adult\_upperExtremities\_Wrist-Supine Feet 1st

Cor GRE T2*	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>	
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>Gradient Echo</i>
	Coil Configuration	<i>GEM Flex Small Array</i>	Imaging Options	<i>NPW</i>
	Plane	<i>OBLIQUE</i>	<b>SCAN RANGE</b>	
	Series Description	<i>Cor GRE T2*</i>	FOV	<i>10.0</i>
	<b>SCAN TIMING</b>		Slice Thickness	<i>2.0</i>
	Flip Angle	<i>33</i>	Slice Spacing	<i>0.5</i>
	TE	<i>Min Full</i>	<b>ACQ TIMING</b>	
	Number of Echoes	<i>1</i>	Freq	<i>256</i>
	TR	<i>700.0</i>	Phase	<i>160</i>
	Receiver Bandwidth	<i>12.50</i>	Freq DIR	<i>S/I</i>
	<b>IMAGE ENHANCE</b>		NEX	<i>2.00</i>
	Filter Choice	<i>None</i>	Auto Shim	<i>Auto</i>
	<b>GATING/TRIGGER</b>		Phase Correction	<i>No</i>
	Auto Trigger Type	<i>Off</i>	<b>FMRI</b>	
	<b>MULTI-PHASE</b>		PSD Trigger	<i>Internal</i>
	Seperate Series	<i>0</i>	View Order	<i>Bottom/Up</i>
	Mask Phase	<i>0</i>	# of Repetitions REST	<i>0</i>
	Mask Pause	<i>0</i>	# of Repetitions ACTIVE	<i>0</i>
	<b>DIFFUSION</b>		<b>SAT</b>	
	Recon All Images	<i>On</i>	Tag Type	<i>None</i>
	<b>CONTRAST</b>		Fat/Water Saturation	<i>Fat</i>
	Contrast Yes/No	<i>No</i>	<b>TRICKS</b>	
			Pause On/Off	<i>On</i>
		Auto Subtract	<i>0</i>	
		Auto SCIC	<i>Off</i>	

Cor GRE T2\*



Protocol: adult\_upperExtremities\_Wrist-Supine Feet 1st

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>FSE-XL</i>
Coil Configuration	<i>GEM Flex Small Array</i>	Imaging Options	<i>NPW, Seq, EDR, TRF, Fast, ZIP512, FR, IrP</i>
Plane	<i>OBLIQUE</i>	<b>SCAN RANGE</b>	
Series Description	<i>Sag STIR</i>	FOV	<i>10.0</i>
<b>SCAN TIMING</b>		Slice Thickness	<i>3.0</i>
Flip Angle	<i>160</i>	Slice Spacing	<i>0.5</i>
TE	<i>34.0</i>	<b>ACQ TIMING</b>	
Number of Echoes	<i>1</i>	Freq	<i>256</i>
TR	<i>3992.0</i>	Phase	<i>160</i>
TI	<i>130</i>	Freq DIR	<i>S/I</i>
Echo Train Length	<i>10</i>	Fat Shift DIR	<i>Normal (I)</i>
Receiver Bandwidth	<i>35.71</i>	NEX	<i>4.00</i>
<b>IMAGE ENHANCE</b>		Auto Shim	<i>Auto</i>
Filter Choice	<i>None</i>	Phase Correction	<i>Yes</i>
<b>GATING/TRIGGER</b>		<b>USER CVS</b>	
Auto Trigger Type	<i>Off</i>	User CV8	<i>1.00</i>
<b>FMRI</b>		User CV22	<i>1.00</i>
PSD Trigger	<i>Internal</i>	TR Min	<i>2500.0</i>
View Order	<i>Bottom/Up</i>	TR Max	<i>8000.0</i>
# of Repetitions REST	<i>0</i>	<b>MULTI-PHASE</b>	
# of Repetitions ACTIVE	<i>0</i>	Seperate Series	<i>0</i>
<b>SAT</b>		Mask Phase	<i>0</i>
Tag Type	<i>None</i>	Mask Pause	<i>0</i>
<b>TRICKS</b>		<b>DIFFUSION</b>	
Pause On/Off	<i>On</i>	Recon All Images	<i>On</i>
Auto Subtract	<i>0</i>	<b>CONTRAST</b>	
Auto SCIC	<i>Off</i>	Contrast Yes/No	<i>No</i>

Sag STIR

Sag STIR

Protocol: adult\_upperExtremities\_Wrist-Supine Feet 1st

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>FSE-XL</i>
Coil Configuration	<i>GEM Flex Small Array</i>	Imaging Options	<i>NPW, EDR, TRF, Fast</i>
Plane	<i>OBLIQUE</i>	<b>SCAN RANGE</b>	
Series Description	<i>Cor T1</i>	FOV	<i>10.0</i>
<b>SCAN TIMING</b>		Slice Thickness	<i>3.0</i>
Flip Angle	<i>120</i>	Slice Spacing	<i>0.0</i>
TE	<i>Min Full</i>	<b>ACQ TIMING</b>	
Number of Echoes	<i>1</i>	Freq	<i>288</i>
TR	<i>581.0</i>	Phase	<i>224</i>
Echo Train Length	<i>3</i>	Freq DIR	<i>S/I</i>
Receiver Bandwidth	<i>35.71</i>	Fat Shift DIR	<i>Normal (S)</i>
<b>IMAGE ENHANCE</b>		NEX	<i>3.00</i>
Filter Choice	<i>None</i>	Auto Shim	<i>Auto</i>
<b>GATING/TRIGGER</b>		Phase Correction	<i>Yes</i>
Auto Trigger Type	<i>Off</i>	<b>USER CVS</b>	
<b>FMRI</b>		User CV8	<i>2.00</i>
PSD Trigger	<i>Internal</i>	User CV9	<i>1.00</i>
View Order	<i>Bottom/Up</i>	User CV22	<i>1.00</i>
# of Repetitions REST	<i>0</i>	TR Min	<i>350.0</i>
# of Repetitions ACTIVE	<i>0</i>	TR Max	<i>800.0</i>
<b>SAT</b>		<b>MULTI-PHASE</b>	
Tag Type	<i>None</i>	Seperate Series	<i>0</i>
<b>TRICKS</b>		Mask Phase	<i>0</i>
Pause On/Off	<i>On</i>	Mask Pause	<i>0</i>
Auto Subtract	<i>0</i>	<b>DIFFUSION</b>	
Auto SCIC	<i>Off</i>	Recon All Images	<i>On</i>
		<b>CONTRAST</b>	
		Contrast Yes/No	<i>No</i>

Cor T1

Cor T1

Protocol: adult\_upperExtremities\_Wrist-Supine Feet 1st

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>Gradient Echo</i>
Coil Configuration	<i>GEM Flex Small Array</i>	Imaging Options	<i>Seq, Fast</i>
Plane	<i>3-PLANE</i>	<b>SCAN RANGE</b>	
Series Description	<i>---Optional---</i>	FOV	<i>20.0</i>
<b>SCAN TIMING</b>		Slice Thickness	<i>5.0</i>
Number of Echoes	<i>1</i>	Slice Spacing	<i>5.0</i>
<b>IMAGE ENHANCE</b>		<b>ACQ TIMING</b>	
Filter Choice	<i>None</i>	Freq	<i>256</i>
<b>GATING/TRIGGER</b>		Phase	<i>128</i>
Auto Trigger Type	<i>Off</i>	Freq DIR	<i>Unswap</i>
<b>MULTI-PHASE</b>		NEX	<i>1.00</i>
Seperate Series	<i>0</i>	Phase FOV	<i>1.00</i>
Mask Phase	<i>0</i>	Auto Shim	<i>Auto</i>
Mask Pause	<i>0</i>	Phase Correction	<i>No</i>
<b>DIFFUSION</b>		<b>FMRI</b>	
Recon All Images	<i>On</i>	PSD Trigger	<i>Internal</i>
<b>CONTRAST</b>		Slice Order	<i>Interleaved</i>
Contrast Yes/No	<i>No</i>	View Order	<i>Bottom/Up</i>
		# of Repetitions REST	<i>0</i>
		# of Repetitions ACTIVE	<i>0</i>
		<b>SAT</b>	
		Tag Type	<i>None</i>
		<b>TRICKS</b>	
		Pause On/Off	<i>On</i>
		Auto Subtract	<i>0</i>
		Auto SCIC	<i>Off</i>

---Optional---

---Optional---

Protocol: adult\_upperExtremities\_Wrist-Supine Feet 1st

Cor T1 PROP	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>	
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>PROPELLER</i>
	Coil Configuration	<i>GEM Flex Small Array</i>	Imaging Options	<i>EDR, TRF, Fast</i>
	Plane	<i>OBLIQUE</i>	<b>SCAN RANGE</b>	
	Series Description	<i>Cor T1 PROP</i>	FOV	<i>14.0</i>
	<b>SCAN TIMING</b>		Slice Thickness	<i>3.0</i>
	Flip Angle	<i>90</i>	Slice Spacing	<i>1.0</i>
	Number of Echoes	<i>1</i>	<b>ACQ TIMING</b>	
	TR	<i>500.0</i>	Freq	<i>288</i>
	Echo Train Length	<i>5</i>	Freq DIR	<i>S/I</i>
	Receiver Bandwidth	<i>50.00</i>	Fat Shift DIR	<i>Normal (I)</i>
	<b>IMAGE ENHANCE</b>		NEX	<i>3.00</i>
	Filter Choice	<i>None</i>	Auto Shim	<i>Auto</i>
	<b>GATING/TRIGGER</b>		Phase Correction	<i>No</i>
	Auto Trigger Type	<i>Off</i>	<b>USER CVS</b>	
	<b>FMRI</b>		User CV5	<i>1.00</i>
	PSD Trigger	<i>Internal</i>	User CV16	<i>1.00</i>
	Slice Order	<i>Interleaved</i>	<b>MULTI-PHASE</b>	
	View Order	<i>Bottom/Up</i>	Seperate Series	<i>0</i>
	# of Repetitions REST	<i>0</i>	Mask Phase	<i>0</i>
	# of Repetitions ACTIVE	<i>0</i>	Mask Pause	<i>0</i>
	<b>SAT</b>		<b>DIFFUSION</b>	
	Tag Type	<i>None</i>	Recon All Images	<i>On</i>
	<b>TRICKS</b>		<b>CONTRAST</b>	
	Pause On/Off	<i>On</i>	Contrast Yes/No	<i>No</i>
	Auto Subtract	<i>0</i>		
	Auto SCIC	<i>Off</i>		

Cor T1 PROP

Protocol: adult\_upperExtremities\_Wrist-Supine Feet 1st

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Feet First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	PROPELLER
Coil Configuration	GEM Flex Small Array	Imaging Options	EDR, TRF, Fast, FR, ARC, IrP
Plane	OBLIQUE		
Series Description	Cor/Ax STIR PROP		
SCAN TIMING		SCAN RANGE	
Flip Angle	160	FOV	12.0
Number of Echoes	1	Slice Thickness	3.0
TR	4337.0	Slice Spacing	0.3
TI	150		
Echo Train Length	19	ACQ TIMING	
Receiver Bandwidth	50.00	Freq	260
		Freq DIR	S/I
		Fat Shift DIR	Normal (S)
		NEX	3.00
		Auto Shim	Auto
		Phase Correction	No
IMAGE ENHANCE		USER CVS	
Filter Choice	None	User CV5	1.00
		TR Min	3500.0
		TR Max	11000.0
GATING/TRIGGER		MULTI-PHASE	
Auto Trigger Type	Off	Seperate Series	0
		Mask Phase	0
		Mask Pause	0
FMRI		DIFFUSION	
PSD Trigger	Internal	Recon All Images	On
View Order	Bottom/Up		
# of Repetitions REST	0	CONTRAST	
# of Repetitions ACTIVE	0	Contrast Yes/No	No
SAT			
SAT Location	I		
Tag Type	None		
TRICKS			
Pause On/Off	On		
Auto Subtract	0		
Auto SCIC	2		

Cor/Ax STIR PROP

Cor/Ax STIR PROP