

PATIENT NAME: [patient name]

EXAM DATE: [exam date] **AGE:** [age] **SEX:** [sex] **I.D.#:** [patient I.D.]

CT CHEST WITH CONT	
CONSENT FORMS	Consent Form Pregnancy Status Form (female patients)
ORAL CONTRAST	NONE
IV CONTRAST	2.5 CC/SEC FOR 75 CC/NS 50 CC
POSITIONING	SUPINE; FEET FIRST; ARMS UP & OUT OF WAY; INSPIRATION BREATH-HOLD CENTER AT STERNAL NOTCH
SCOUT	S 60/I 300, 120 KVP @ 10 MA, AP & LAT
TECHNIQUE - 16 SLICE	HELICAL; 2.5 MM, STD W400/L40 KVP 120 @ AUTO/SMART MA - 40% ASIR SFOV - LGE BODY/DFOV - FIT TO ANATOMY ROI AORTA AT THE LEVEL OF HEART 90 HU WITH MONITOR DELAY OF 10 SEC SCAN FROM APEX THROUGH ADRENALS
TECHNIQUE - 64 SLICE	HELICAL; 2.5 MM, STD W400/L40 KVP 120 @ AUTO/SMART MA - 40% ASIR SFOV - LGE BODY/DFOV - FIT TO ANATOMY ROI AORTA AT LEVEL OF THE HEART SCAN FROM APEX THROUGH ADRENALS
TECHNIQUE - 128 SLICE	HELICAL; 2.5MM, STD W400/L40 KVP 120 @ AUTO/SMART MA - 40% ASIR SFOV - LGE BODY/DFOV - FIT TO ANATOMY DELAY 40 SEC SCAN FROM APEX THROUGH ADRENALS
TECHNIQUE - SIEMENS	SPIRAL; 3 MM /19.2 / 1:0.8, STD, 0.5 SEC KVP 100, CAREdose 4D, CARE KVP off, AUTO MA, Dose saving optimized 7, Lower limit max mAs 25% DELAY 40 SEC, SCAN FROM APEX THROUGH ADRENALS
RECONS	2.5 MM, LUNG W1400/L-600 1.25 MM, STD W400/L40 0.625 MM, STD W400/L40 (DMPR) 16 / 128 SLICE: 1.25 MM/INT 0.625 MM, STD W400/L40

	Siemens- 0.625 MM STD W 400/L40 3 MM LUNG W 1400/-600
REFORMATIONS	CORONAL & SAGITTAL AVG 2.5MM/2.5MM IN STD W400/L40 AXIAL THICK MIP 6MM X 3MM, STD W1400/L-600 Siemens- Coronal 3 MM, Lung W1400/L-600
PACS	SCOUT 2.5 MM - LUNG 2.5 MM - CHEST CHEST COR & CHEST SAG AX LUNG MIP CHEST 1.25 Siemens- 3 MM Lung & Chest w
CHARGE	CT CHEST WITH CONTRAST USED
REMARKS	