

<b>CT KNEE POST ARTHROGRAM</b>	Indication: POST ARTHROGRAM R/O TEAR
FORMS	Consent Form
ORAL CONTRAST	NONE
IV CONTRAST	NO IV CONTRAST. CONTRAST IS INJECTED INTO THE JOINT IN XRAY.
POSITIONING	SUPINE, FEET FIRST, AFFECTED KNEE IN CENTER OF TABLE; WITH TOES STRAIGHT UP CENTER MID KNEE
SCOUT	S150/I 150, 120 KVP @ 10 MA, AP S150 / I 150, 120 KVP @ 10 MA, LAT
TECHNIQUE - 16 SLICE	HELICAL; 0.625 MM/INT 0.25 MM, STD W400/L40 120 KVP @ AUTO MA -30% ASIR DFOV- 18 CM (FIT TO ANATOMY)  IMMOBILIZE KNEE, SCAN ENTIRE KNEE INCLUDE ALL OF KNEE JOINT AND ANY CONTRAST OUTSIDE OF THE JOINT CAPSULE.
TECHNIQUE - 64 SLICE	HELICAL; 0.625 MM, 120 KVP @ 335 MA 30% ASIR STD W400/L40 DFOV-18 CM (FIT TO ANATOMY) IMMOBILIZE KNEE; SCAN ENTIRE KNEE INCLUDE ALL OF KNEE JOINT AND ANY CONTRAST OUTSIDE OF THE JOINT CAPSULE.
TECHNIQUE – 128 SLICE	
RECONS	1.25 MM BONE PLUS W2500/L250 1.25 MM STD W400/L40 0.625 MM, STD W/400/L40 (DMPR)
REFORMATIONS	OBLIQUE CORONAL & OBLIQUE SAGITTAL TO THE FEMORAL CONDYLES - AVG 0.6 MM THICK/1.5 MM SPACING, W400/L40, BONE PLUS & STD
PACS	SCOUT 1.25 MM – KNEE POST ARTH 1.25 MM - STD EXT COR & EXT SAG – BONE & STD
CHARGE	KNEE POST ATHROGRAM
REMARKS	