

CT MAKO HIP	MAKO THA PROTOCOL
CONSENT FORMS	Pregnancy Status Form (Female Patients)
ORAL CONTRAST	NO
IV CONTRAST	NO
POSITIONING	SUPINE; FEET FIRST, ARMS OVER CHEST; TOES STRAIGHT UP, ENSURE PT IS IN TRUE SUPINE POSITION BY PALPATING ANTERIOR SUPERIOR ILIAC SPINES AND COMPARING RELATIVE HEIGHT ABOVE THE CT TABLE. ALIGN LONGITUDINAL AXIS OF BODY WITH LONGITUDINAL AXIS OF CT TABLE. REMOVE ALL CLOTHING WITH METAL. MAKE SURE LEGS ARE IN THE SAME POSITION. ***NO METAL BAR USED FOR THIS STUDY
SCOUT	S50/I1000, 120 KVP @ 10 MA, AP S50/I1000 120 KVP @ 10 MA, LAT SCOUTS ARE USED FOR MEASUREMENT
TECHNIQUE- ALL	HELICAL; W2500/L250 120-140 KVP @ 250 MAX AUTO/SMART MA, BONE, 30% ASIR NO BREATH HOLD SFOV - LGE BODY / DFOV- FIT TO INCLUDE PELVIS & BOTH KNEES SCAN PELVIS: 0.625MM, ENTIRE PELVIS THRU TO MID THIGH SCAN BOTH KNEES: 5 MM, JOINT LINES BETWEEN FEMUR AND TIBIA INCLUDING 10 CM PROXIMAL TO JOINT LINE OF FEMUR
RECONS	NONE
REFORMATIONS	NONE
PACS	SCOUT MAKO HIP
CHARGE	CT MAKO HIP RIGHT OR LEFT
REMARKS	SEND TO MAKO DESTINATION