

CT MAKO KNEE	MAKO PKA PROTOCOL
CONSENT FORMS	Pregnancy Status Form (Female Patients)
ORAL CONTRAST	NO
IV CONTRAST	NO
POSITIONING	<p>SUPINE; FEET FIRST, ARMS OVER CHEST; TOES STRAIGHT UP WITH METAL ROD TUCKED NEXT TO HIP DOWN TO ANKLE. ATTACH VELCRO ON THE METAL ROD AROUND THIGH AND ANKLE.</p> <p>*SMALL TOWEL UNDER KNEE.</p> <p>REMOVE ANY CLOTHING WITH METAL.</p> <p>** IF PT HAS METAL IN THE OPPOSITE KNEE, MAKE SURE AND SEPARATE KNEES TO DECREASE ARTIFACT FROM THE METAL.</p>
SCOUT	<p>S50/I1500, 120 KVP @ 10 MA, AP</p> <p>S50/I1500 120 KVP @ 10 MA, LAT</p> <p>SCOUTS ARE USED FOR MEASUREMENT</p>
TECHNIQUE- ONE SERIES W/ THREE GROUPS	<p>HELICAL; W2500/L250</p> <p>120 KVP @ AUTO/SMART MA, BONE, 30% ASIR APPLIED; NO BREATH HOLD</p> <p>SFOV - LGE BODY / DFOV –</p> <p> HIP- DON'T EXCEED 50CM</p> <p> KNEE – DON'T EXCEED 25CM</p> <p> ANKLE- DON'T EXCEED 50CM</p> <p>**INCLUDING METAL ROD IN EVERY IMAGE</p> <p>SCAN HIP: 5 MM, FROM GREATER SCIATIC NOTCH THRU TO LESSER TROCHANTER</p> <p>SCAN KNEE: 0.625MM, 10 CM ABOVE & 10 CM BELOW THE KNEE.</p> <p>SCAN ANKLE: 5 MM, THROUGH ANKLE JOINT</p>
RECONS	NONE
REFORMATIONS	NONE
PACS	SCOUT MAKO KNEE
CHARGE	CT MAKO KNEE RIGHT OR LEFT
REMARKS	SEND TO MAKO DESTINATION