

CT FACIAL BONES WITHOUT PEDIATRIC	Assess fractures of face and orbits
CONSENT FORMS	Consent Form
ORAL CONTRAST	NONE
IV CONTRAST	NONE UNLESS PHYSICIAN ORDERS WITH CONTRAST
POSITIONING	SUPINE; HEAD FIRST IN HEAD HOLDER; ARMS DOWN BY SIDE; CENTER PERPENDICULAR TO OM BASELINE
SCOUT	S 150/I 150, 100-120 KVP @ 10 MA, AP & LAT
TECHNIQUE - 16 SLICE	HELICAL, 5.0 MM, BONE W3000/L300, PITCH 0.983 SEE TECHNIQUE CHARTS BELOW SCAN FROM SUPRA-ORBITAL RIM DOWN THROUGH MANDIBLE
TECHNIQUE - 64 SLICE	HELICAL; 0.625 MM, W3000/L300, PITCH 0.969 SEE TECHNIQUE CHARTS BELOW SCAN FROM SUPRA-ORBITAL RIM DOWN THROUGH MANDIBLE
TECHNIQUE – 128 SLICE	
TECHNIQUE - SIEMENS	
RECONS	16 SLICE-1.25 MM STD 16 SLICE-1.25/INT 0.625 MM BONE W3000/L300 -DMPR 64 SLICE-1.25 MM STD 16 & 64 SLICE-1.25 MM BONE W3000/L300
REFORMATIONS	SAGITTAL & CORONAL AVG 0.6 MM THICK / 1.5 MM SPACING, W3000/L300
PACS	SCOUT 1.25 MM - FACIAL WO 1.25 MM - STD SINUS COR & SINUS SAG
CHARGE	FACIAL BONES WITHOUT
REMARKS	

