

CT SINUS LANDMARX	MEDTRONIC NAVIGATION, FUSION
CONSENT FORMS	Consent to Contrast Material
ORAL CONTRAST	NONE
IV CONTRAST	NONE
POSITIONING	SUPINE; HEAD FIRST IN ON A BLANKET/FOAM HEADREST ARMS DOWN BY SIDE; CENTER PERPENDICULAR TO INFRAORBITAL BASELINE;
SCOUT	S150 / I 150, 120 KVP @ 10 MA, AP & LAT
TECHNIQUE - 64 SLICE	HELICAL; 1.25 MM, W400/L40 140 KVP @ 115 MA, 30% ASIR FOV-25 CM SCAN FROM SUPERIOR ASPECT OF HORIZONTAL PORTION OF THE MANDIBLE THROUGH VERTEX. INCLUDE EARS, MAXILLARY TEETH, TIP OF NOSE AND POSTERIOR SKULL.
RECONS	64 SLICE: 0.625 MM, BONE W3000/L300 (DMPR) 1.25 MM BONE W3000/L300
REFORMATIONS	CORONAL AND SAGITTAL - AVG 0.6 MM THICK/1.5 SPACING, W3000/L300
PACS	SCOUT 1.25 MM – LANDMARX SINUS 1.25 MM - BONE SINUS COR & SINUS SAG
CHARGE	CT LANDMARX SINUS
REMARKS	OHSU DOCTORS- PUSH LANDMARX STD SERIES OVER TO OHSU DOCTORS. LOCAL DOCTORS- WILL CALL MEDICAL RECORDS FOR A DISC TO BE MADE OF LANDMARX STD SERIES.