

<b>PEDIATRIC SINUS</b>	
CONSENT FORMS	Consent Form
ORAL CONTRAST	NONE
IV CONTRAST	NONE
POSITIONING	SUPINE; HEAD FIRST IN HEAD REST; RAISE BODY LEVEL WITH HEAD REST; USE VELCRO STRAPS AND IMMOBILIZING DEVICES; CENTER AT OM BASELINE, PERPENDICULAR TABLE
SCOUT	S 120 / I 100, 120 KVP@ 10 MA, LAT & AP
TECHNIQUE	HELICAL; 5 MM / 9.37/0.938;1, BONE W3000/L300 120 KVP @ 200 MA, NI 6.0 - 30% ASIR SFOV - HEAD / DFOV - 22 CM SCAN FROM TOP OF FRONTAL SINUS THROUGH CANTHUS
RECONS	1.25 MM, BONE PLUS W3000/L300 2.5 MM, STD W400/L40 0.625 MM, BONE W3000/L300
REFORMATIONS	CORONAL & SAGITTAL - AVG 0.6 MM THICK/ 1.5 MM SPACING, W3000/L300
PACS	SCOUT 1.25 MM - SINUS 2.5 MM - STD SINUS COR & SINUS SAG
CHARGE	SINUS COMPLETE
REMARKS	WORK WITH PARENTS, IF PRESENT, TO KEEP CHILD CALM AND TO HELP IMMOBILIZE HEAD.